

Therapy optimization on the basis of regular FeNO measurement

Optimization of therapy based on regular FeNO home measurement by quantifying Type 2 inflammation

Dr. med. Hartmut Timmermann, Pneumologikum Hamburg, Valentinskamp 24, 20354 Hamburg

Summary

The case study refers to a participant of the ongoing FeNO@home study. The patient with asthma known since childhood had reported good asthma control in the last few years. His lung function was normal. He performed FeNO home measurements over a period of 12 weeks under his usual medication. We found unexpectedly high FeNO values confirming an ongoing Type 2 inflammation. He was instructed to double his ICS dose with normalizing of FeNO. **The case illustrates how frequent FeNO home measurements can lead to optimized asthma treatment.**

Case history

The case described in this report relates to a 44-year-old male patient who was diagnosed with asthma in childhood. He is a non-smoker. He has allergies to grass and pollen, house dust mites, and cat epithelium. The patient presents a known comorbidity of allergic rhinitis.

He showed controlled asthma symptoms according to the Asthma Control Test (ACT score 21 units). The lung function was considered as normal (FEV1 of 104%). Measurements of Fractional exhaled Nitric Oxide (FeNO) in medical practice indicated a high FeNO-level (95 ppb).

The patient received GINA (Global Initiative of Asthma) treatment step 2-3 consisting of low-dose Symbicort (2x1).

Problem statement

What is the best approach for patients with well controlled asthma and high FeNO values confirming an ongoing Type 2 inflammation?

Investigation

The patient was included in the FeNO@home study¹ outside pollen season. The aim of the study was to investigate whether regular FeNO home measurements had an impact on patient compliance or behavior, variability of FeNO values over a longer period, correlation of FeNO values with symptoms, identification of asthma triggers, and treatment decisions. In this multicenter study, adult patients with diagnosed asthma performed FeNO measurements over a period of 12 weeks using the Vivatmo *me* measurement device for home use. They continued to take their currently prescribed asthma treatment, which could also be adapted. Daily symptoms, use of asthma medication, potential exacerbations, and Peak Expiratory Flow (PEF) were recorded in the device-associated Vivatmo *app*. After 12 weeks, the study ended with a final assessment of asthma control, symptoms, and lung function.

Results and treatment

The medication was adjusted during the study. With normal lung function and a high ACT score, we found unexpectedly high FeNO values confirming an existing Type 2 inflammation. During the first month, the FeNO values were in a very high range (50-80 ppb), leading to an increase of the dose (2x2). After increasing the ICS dose, the majority of FeNO levels were below the 25-ppb-threshold from the ATS guideline². The low FeNO values were accompanied by less coughing and an improvement in general health of the patient. The final assessment in the medical practice indicated normal lung function (FEV1 4,72 l, 103% of expected; FVC 5,49 l; FeNO 18 ppb), and an ACT score of 23 units.

- 4 Guideline for the Diagnosis and Treatment of Asthma – Addendum 2020 Guideline of the German Respiratory Society and the German Atemwegsliga in Cooperation with the Paediatric Respiratory Society and the Austrian Society of Pneumology. *Pneumologie* 2021;75(03):191-200
- 5 Beasley R, Holliday M, Reddel HK, et al. Controlled trial of budesonide-formoterol as needed for mild asthma. *N Engl J Med* 2019; 380: 2020-2030
- 6 Hardy J, Baggott C, Fingleton J, et al. Budesonide-formoterol reliever therapy versus maintenance budesonideplus terbutaline reliever therapy in adults with mild to moderate asthma (PRACTICAL): a 52-week, open-label, multicentre, superiority, randomised controlled trial. *Lancet* 2019; 394: 919-928.
- 7 Diagnosis and treatment of asthma: a guideline for respiratory specialists 2023 published by the German Respiratory Society (DGP). *Pneumologie* 2023. DOI 10.1055/a-2070-2135
- 8 Pavord ID, Beasley R, Agusti A, et al. After asthma: redefining airways diseases. *Lancet* 2018. DOI: 10.1016/S0140-6736(17)30879-6.